


KIRIYAMA Sublimation Apparatus for Very Small Quantity



www.kiriyama.co.jp

W77A-3 type sublimation apparatus for a very small quantity

We developed the convenient sublimation apparatus specialized for a very small quantity base on the Mill-type device. The design is brushed up as the best.

<Features>

By combining the joint portion and the center condensing pipe, it is very simple and easy to wash and the endurance is enhanced.



Compact condensing unit and sample tube



Condensing unit with no dead space



Example of experiment



**Purification by sublimation
(crystal like needle)**



Left W77A-3-2

Right W77A-3-1

We are preparing two types
W77A-3-1(joint size 19/38)19,000yen
W77A-3-2(joint size 24/40)22,000yen
which has larger caliber.

※ Please be aware that line-up of products and shapes may change in the future.

KIRIYAMA GLASS WORKS COMPANY

2-31-11, Higashi-nippori, Arakawa-ku, Tokyo, Japan, 116-0014

URL <http://kiriyaama.co.jp> TEL 03(3802)0005

e-mail : info@kiriyaama.co.jp FAX 03(3801)1170